## IS-III active



2019.01 E-Ver.03



#### Contents

IS-III active Benefits	P 06
IS-III active Features	P 08
IS-III active Line Up	P11
Clinical Cases of IS-III active	P 12
Neo Surgical Kit	P <b>14</b>
Opening the implant ampoule	P 19
IS-III active Drilling Protocol	P 20



## $\mathsf{V} \mathsf{H} \mathsf{Y} \mathsf{IS-III} \mathsf{active}^?$

IS-III active implant is structured to maximize initial stability and facilitate faster osseointegration with its scientifically proven SLA surface and fixture body design.





Ð

### 01 Connection

- A. Thicker Platform
- **B**. Anti-screw Loosening
- C. Abutment Compatibility

### )2 Design

- D. Platform Microgroove
- E. Magic Threads
- F. Deep & Wide Pitch
- **G**. Cutting Edge



- H. S.L.A. Surface
- I. Cell Adhesion Ability





### IS-III active Benefits



### 01 Connection

Anti-screw Loosening ► Two Connection Points ► Eliminate screw fracture

Abutment Compatibility ► Compatible with IS type ► Conical 11° / Internal 2.5 Hex

### 02 Design

Platform Microgroove 🕨 Enhanced Soft Tissue Sealing

Minimize bone loss

Deep&Wide Pitch ► Reduced Bone Compression

Optimal for Osseointegration

Wider Cutting Edge ► Improved Self-tapping Ability

Maximize initial stability

Magic Threads ► Endure Vertical/Lateral Force ► Maximize initial stability

### )3 Surface

Improved Surface > Increased Surface Area

► Facilitate faster osseointegration

Greater Cell Adhesion Ability ► More Cell adhesion

- ► Facilitate faster osseointegration
- - Predictable Implant Placement
  - Successful Primary & Secondary Stability
  - Faster Patient Recovery & Masticatory Function

# IS-III active Features Platform & Connection





Minimize Microgroove design at the upper platform of the fixture enhances soft tissue sealing, thus prevents bone loss.



Platform microgroove



Enhanced soft tissue barrier seal



Minimize bone loss through soft tissue integration and optimized soft tissue seal

The coronal area of the fixture is also S.L.A. surface treated and takes a bevel border with open threaded design. These features facilitate osseointegration to crestal bone level, as well as minimize bone loss and maintain bone level.





Open threaded bevel coronal

Minimize bone loss & maintain bone level



Successful osseointegration to bone level

Thicker connection through Increased platform thickness.



Maintains connection thickness over 3mm



Increased strength of connection

#### Stronger Connection



### S.L.A. Surface



The new S.L.A. Surface with 40% greater surface area and 50% more cell adhesion promotes faster osseointegration.





Improved processing technique of the S.L.A. Surface



40 percent increase in surface area



Reduced osseointegration time (50 percent increase in cell adhesion)

### Wide Cutting Edge





Wider cutting edge and enlarged surface area enhances initial fixation and offers clinicians more stable implant placement.



Doubled cutting edge surface



Improved Self-tapping ability while minimizing bone compression



Maximized initial fixation (AnyTime Loading)

# IS-III active Features Deep&Wide Pitch



### Surgical Kit

Optimum Pitch for Osseointegration.



Increase in thread pitch to 0.9



Minimal bone compression (Prevent bone necrosis)



Provide optimal condition for osseointegration

More accessibility with improved cutting force of the surgical drills, now available in two different lengths.





Clinicians decide the loading time by utilizing either the cortical drill or the cortical tap according to the patient's bone density and oral conditions.

#### Cortical Drill

Cortical Tap

Utilized for Delayed Loading by drilling the crestal cortical bone.

Utilized for Immediate (Any-Time) Loading by tapping

the crestal cortical bone.



Delayed Loading



Immediate Loading (AnyTime Loading)



### 10/11

### IS-III active Line Up

Platform

Diameter



\* Coverscrew not included.

### Clinical Cases of IS-III active

Case 1



Pre-op panorama (#46, 47)



Intra-oral photograph



Flap reflection



Bone trimming for osteotomy



Drilling & Cortical tapping



After cortical tapping



IS-III active placement in #46 and #47



ITV of 40Ncm for both sites



Healing abutment & suture



Post-op panorama on the day of surgery



Final restorations after 5 months



6-months follow-up radiograph

#### Case 2



Pre-op panorama (#36)



Intra-oral photograph



Flap reflection



Drilling & Cortical Tapping



Removing fixture from the ampoule



IS-III active placement in #36



ITV of 40Ncm



Healing abutment & Suture



Post-op panorama (#36)



Final restoration after 2 months



7-months follow-up radiograph

### Neo Surgical Kit



14/15

#### Point Lindemann Drill





39.4	33.4	F 02.2	1025					
1.2			Ņ	7.3	8.5	10	11.5	13

Initial Drill

Diameter	Туре	Product name
Ø2.2	Short	TSD22CS
Ø2.2	Long	TSD22CL

Twist Drill



Diameter	Туре	Product name
Ø3.0	Short	TSD30CS
Ø3.0	Long	TSD30CL
Ø3.5	Short	TSD35CS
Ø3.5	Long	TSD35CL
Ø4.0	Short	TSD40CS
Ø4.0	Long	TSD40CL
Ø4.5	Short	TSD45CS
Ø4.5	Long	TSD45CL

Stopper

Drilling Length

	3.0	4.0	5,0	6,0	0.0	7.3	8.5	40.0	411.5	43.0
Stopper	3.0	4.0	5.0	6.0	6.6	7.3	8.5	10.0	11.5	13
Drilling Length(mm)	4.2	5.2	6.2	7.2	7.8	8.5	9.7	11.2	12.7	14.2
Product name	DS030C	DS040C	DS050C	DS060C	DS066C	DS070C	DS085C	DS100C	DS115C	DS130C

Ø

### Neo Surgical Kit



Diameter	Product name
Ø3.65	ISCD35F
Ø4.2	ISCD40F
Ø4.4	ISCD45F
Ø4.9	ISCD50F

#### **IS Fixture Driver**



Length	Product name
Ratchet (Short)	ISFD10R
Ratchet (Long)	ISFD15R
Contra Angle (Short)	ISFD05C
Contra Angle (Long)	ISFD05CL

#### Cortical Tap



Diameter	Product name
Ø3.5	ISTD38S
Ø4.0	ISTD43S
Ø4.5	ISTD45S
Ø5.0	ISTD50S

#### Connector



Length	Product name
Short	RC10
Long	RC15

Product name CAA00

#### **Direction Pin**



Diameter	Product name
Ø3.5	DPIS35C
Ø4.5	DPIS45C

#### **Parallel Pin**



Length	Product name
7.0mm	PP07F
8.5mm	PP08F
10.0mm	PP10F

### 16/17



### Neo Surgical Kit



Twist Drill

Twist Drill + Contra Angle Twist Drill + Extention + Contra Angle



### Opening the implant ampoule



Remove the square-shaped ampoule from the blister.



Turn the lid to open the ampoule.



Remove the inner circular ampoule from the outer square-shaped ampoule.



Drop the inner ampoule onto the operating table.



Remove the safety cap (A cover screw can be found inside the cap).



Hold the sides of the ampoule when removing the cap. Must be cautious not to grip on the clip. (Opening of the clip will cause the fixture to fall into the ampoule.)



Hold the upper part of the clip and connect the fixture driver to the implant.



Simultaneously, push the lower part of the clip for clip opening and lift the implant out of the ampoule.

### IS-III active Drilling Protocol

#### IS-III active Fixture Ø3.5 X 10mm (D1/D2 bone)



In soft(D4) bone, use Ø 2.2 initial drill as the final drill



#### IS-III active Fixture Ø4.0 X 10mm (D1/D2 bone)

In soft(D4) bone or in condition of getting initial fixation at implant apex, Ø3.0 twist drill is the final drill

#### **Drilling Speed & Torque**

Point Lindemann, Initial Drill, Twist Drill : 1,200rpm / 35~45Ncm Cortical Tap : 50rpm / 50Ncm Cortical Drill : 1200rpm / 50Ncm (Conventional Loading case)

#### IS-III active Fixture Ø4.5 X 10mm (D1/D2 bone)



In soft(D4) bone or in condition of getting initial fixation at implant apex, Ø3.5 twist drill is the final drill



#### IS-III active Fixture Ø5.0 X 10mm (D1/D2 bone)

In soft(D4) bone or in condition of getting initial fixation at implant apex, Ø4.0 twist drill is the final drill





#### History of Neobiotech

Mar. 2017	Ridge Wider Kit
Feb. 2017	T-brush
Sep.2016	IS-III active
Jul. 2016	EZ GBR System
May 2015	Encoded Healing abutment
Apr. 2015	CAMeleon cs
May 2014	World Class 300
Dec. 2013	Manufactured CAMeleon
Nov. 2013	EB-II active
Oct. 2013	SinusAll Kit
	PickCap Impression Kit
Jun. 2013	IT-II active
Oct. 2012	Prosthetic Kit / Accessory Kit
Jun. 2012	Neoguide system
Mar. 2012	GBR Kit
Oct. 2011	IS-II active, Quicktight
Jun. 2011	IS-II, S-mini & ACM
Oct. 2010	CTi – mem
Feb. 2010	SR Kit
Jun. 2009	FR Kit
Mar. 2009	Wide Implant
Nov. 2008	CMI IS implant
Jul. 2008	SLA-Kit
Mar. 2008	SCA-Kit
Mar. 2008	Obtain the patent of CMI Implant
Sep. 2007	Merged with "Osscare.Co.Ltd"
Jun. 2007	CMI implant(External Type)
Feb. 2007	Change of Management
Jul. 2000	Foundation of "Neobiotech.Co,Ltd,."



Neobiotech Co.,Ltd 10F, 36, Digitalro-27Gil, Guro-Gu, Seoul, Korea, 08381 Tel 82-2-582-2885, Fax 82-2-582-2883 www.neobiotech.com / Made in Korea